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- O'Brien E, Atkins N, Stergiou G, Karpettas N, Parati G, Asmar R, Imai Y, Wang J, Mengden T, Shennan A; Working Group on Blood Pressure Monitoring of the European Society of Hypertension. European Society of Hypertension International Protocol revision 2010 for the validation of blood pressure measuring devices in adults. Blood Press Monit. 2010 Feb;15(1):23-38. Erratum in: Blood Press Monit. 2010 Jun;15(3):171-2.
- O'Brien E, Petrie J, Littler W, de Swiet M, Padfield PL, Altman DG, Bland M, Coats A, Atkins N. An outline of the revised British Hypertension Society protocol for the evaluation of blood pressure measuring devices. J Hypertens. 1993 Jun;11(6):677-9.
- Stergiou GS, Tzamouranis D, Protogerou A, Nasothimiou E, Kapralos C. Validation of the Microlife Watch BP Office professional device for office blood pressure measurement according to the International protocol. Blood Press Monit. 2008 Oct;13(5):299-303.
- Lohmann FW, Eckert S, Verberk WJ. Interarm differences in blood pressure should be determined by measuring both arms simultaneously with an automatic oscillometric device. Blood Press Monit. 2011 Feb;16(1):37-42.
- Saladini F, Benetti E, Masiero S, Palatini P. Accuracy of Microlife WatchBP Office ABI monitor assessed according to the 2002 European Society of Hypertension protocol and the British Hypertension Society protocol. Blood Press Monit. 2011 Oct;16(5):258-61.
- Stergiou GS, Tzamouranis D, Nasothimiou EG, Karpettas N, Protogerou A. Are there really differences between home and daytime ambulatory blood pressure? Comparison using a novel dual-mode ambulatory and home monitor. J Hum Hypertens. 2010 Mar;24(3):207-12.
- Ragazzo F, Saladini F, Palatini P. Validation of the Microlife WatchBP O3 device for clinic, home, and ambulatory blood pressure measurement, according to the International Protocol. Blood Press Monit. 2010 Feb;15(1):59-62.
- Stergiou GS, Giovas PP, Gkinos CP, Patouras JD. Validation of the Microlife WatchBP Home device for self home blood pressure measurement according to the International Protocol. Blood Press Monit. 2007 Jun;12(3):185-8.
- Stergiou GS, Giovas PP, Neofytou MS, Adamopoulos DN. Validation of the Microlife BPA100 Plus device for self-home blood pressure measurement according to the International Protocol. Blood Press Monit. 2006 Jun;11(3):157-60.
- Palatini P, Dorigatti F, Bonso E, Ragazzo F. Validation of Microlife BP W100 wrist device assessed according to the European Society of Hypertension and the British Hypertension Society protocols. Blood Press Monit. 2009 Feb;14(1):41-4.
- Cuckson AC, Reinders A, Shabeeh H, Shennan AH; British Hypertension Society. Validation of the Microlife BP 3BTO-A oscillometric blood pressure monitoring device according to a modified British Hypertension Society protocol. Blood Press Monit. 2002 Dec;7(6):319-24.
- Chung Y, de Greeff A, Shennan A. Validation and compliance of a home monitoring device in pregnancy: microlife WatchBP home. Hypertens Pregnancy. 2009;28(3):348-59.
- Reinders A, Cuckson AC, Lee JT, Shennan AH. An accurate automated blood pressure device for use in pregnancy and pre-eclampsia: the Microlife 3BTO-A. BJOG. 2005 Jul;112(7):915-20.

- 14 Wilton A, De Greef A, Shennan A. Rapid assessment of blood pressure in the obstetric day unit using Microlife MaM technology. *Hypertens Pregnancy*. 2007;26(1):31-7.
- 15 Simões GC, Moreira SR, Kushnick MR, Simões HG, Campbell CS. Postresistance exercise blood pressure reduction is influenced by exercise intensity in type-2 diabetic and nondiabetic individuals. *J Strength Cond Res*. 2010 May;24(5):1277-84.
- 16 Motta DF, Lima LC, Arsa G, Russo PS, Sales MM, Moreira SR, Moraes PK, Almeida WS, Araujo RC, Moraes MR, Pesquero JL, Simões HG, Campbell CS. Effect of type 2 diabetes on plasma kallikrein activity after physical exercise and its relationship to post-exercise hypotension. *Diabetes Metab*. 2010 Nov;36(5):363-8.
- 17 Thompson AM, Eguchi K, Reznik ME, Shah SS, Pickering TG. Validation of an oscillometric home blood pressure monitor in an end-stage renal disease population and the effect of arterial stiffness on its accuracy. *Blood Press Monit*. 2007 Aug;12(4):227-32.
- 18 Masiero S, Saladini F, Benetti E, Palatini P. Accuracy of the Microlife large-extra large-sized cuff (32-52 cm) coupled to an automatic oscillometric device. *Blood Press Monit*. 2011 Apr;16(2):99-102.
- 19 Wiesel J, Fitzig L, Herschman Y, Messineo FC. Detection of atrial fibrillation using a modified microlife blood pressure monitor. *Am J Hypertens*. 2009 Aug;22(8):848-52.
- 20 G. S. Stergiou, A. Destounis, A. Kollias, D. Tzamouranis, N. Karpettas, P. Kalogeropoulos, E. Andreadis- Accuracy of automated oscillometric blood pressure measurement in patients with atrial fibrillation. *J Hypertens* 2011;29(Suppl.A):e2.
- 21 Parati G, Omboni S, Palatini P, Rizzoni D, Bilo G, Valentini M, Rosei EA, Mancia G. Italian society of hypertension guidelines for conventional and automated blood pressure measurement in the office, at home and over 24 hours. *High Blood Press Cardiovasc Prev*. 2008 Oct;15(4):283-310.
- 22 O'Brien E, Asmar R, Beilin L, Imai Y, Mancia G, Mengden T, Myers M, Padfield P, Palatini P, Parati G, Pickering T, Redon J, Staessen J, Stergiou G, Verdecchia P; European Society of Hypertension Working Group on Blood Pressure Monitoring. Practice guidelines of the European Society of Hypertension for clinic, ambulatory and self blood pressure measurement. *J Hypertens*. 2005 Apr;23(4):697-701.
- 23 Pickering TG, Hall JE, Appel LJ, Falkner BE, Graves J, Hill MN, Jones DW, Kurtz T, Sheps SG, Roccella EJ. Recommendations for blood pressure measurement in humans and experimental animals: part 1: blood pressure measurement in humans: a statement for professionals from the Subcommittee of Professional and Public Education of the American Heart Association Council on High Blood Pressure Research. *Circulation*. 2005 Feb 8;111(5):697-716.
- 24 Yarows SA, Patel K, Brook R. Rapid oscillometric blood pressure measurement compared to conventional oscillometric measurement. *Blood Press Monit*. 2001 Jun;6(3):145-7.
- 25 Dourmap-Collas C, Villeneuve F, Jaboureck O, Chantrel F, Hanon O, Girerd X; Club des jeunes hypertensiologues. [Assessment of home blood pressure with a monitor including MAM technology: comparison with the standard monitor]. *Arch Mal Coeur Vaiss*. 2006 Jul-Aug;99(7-8):754-7.
- 26 Palatini P, Dorigatti F, Bonso E, Ragazzo F. Validation of the Microlife BP W200-1 wrist device for blood pressure measurement. *Blood Press Monit*. 2008 Oct;13(5):295-8.
- 27 European Heart Rhythm Association; European Association for Cardio-Thoracic Surgery, Camm AJ, Kirchhof P, Lip GY, Schotten U, Savelieva I, Ernst S, Van Gelder IC, Al-Attar N, Hindricks G, Prendergast B, Heidbuchel H, Alfieri O, Angelini A, Atar D, Colonna P, De Caterina R, De Sutter J, Goette A, Gorenek B, Haldal M, Hohloser SH, Kolh P, Le Heuzey JY, Ponikowski P, Rutten FH. Guidelines for the management of atrial fibrillation: the Task Force for the Management of Atrial Fibrillation of the European Society of Cardiology (ESC). *Eur Heart J*. 2010 Oct;31(19):2369-429. Erratum in: *Eur Heart J*. 2011 May;32(9):1172.
- 28 Zoni-Berisso M, Filippi A, Landolina M, Brignoli O, D'Ambrosio G, Maglia G, Grimaldi M, Ermini G. Frequency, patient characteristics, treatment strategies, and resource usage of atrial fibrillation (from the Italian Survey of Atrial Fibrillation Management [ISAF] study). *Am J Cardiol*. 2013 Mar 1;111(5):705-11.
- 29 Manolis AJ, Rosei EA, Coca A, Cifkova R, Erdine SE, Kjeldsen S, Lip GY, Narkiewicz K, Parati G, Redon J, Schmieder R, Tsioufis C, Mancia G. Hypertension and atrial fibrillation: diagnostic approach, prevention and treatment. Position paper of the Working Group 'Hypertension Arrhythmias and Thrombosis' of the European Society of Hypertension. *J Hypertens*. 2012 Feb;30(2):239-52.
- 30 NICE, National Institute for Health and Care Excellence. Diagnosis and assessment of hypertension. 2013. <http://pathways.nice.org.uk/pathways/hypertension>, 2013
- 31 NICE, National Institute for Health and Care Excellence. WatchBP Home A for opportunistically detecting atrial fibrillation during diagnosis and monitoring of hypertension. 2013. <http://guidance.nice.org.uk/MTG13>
- 32 Battigelli D, Brignoli O, Ermini G, Filippi A, Guilaro B, Giustini SE. Serie Editoriale Disease Management. Società Italiana di Medicina Generale (SIMG). Fibrillazione Atriale in Medicina Generale. 2013.
- 33 Wiesel J, Wiesel D, Suri R, Messineo FC. The use of a modified sphygmomanometer to detect atrial fibrillation in outpatients. *Pacing Clin Electrophysiol*. 2004 May;27(5):639-43.
- 34 Wiesel J, Fitzig L, Herschman Y, Messineo FC. Detection of atrial fibrillation using a modified microlife blood pressure monitor. *Am J Hypertens*. 2009 Aug;22(8):848-52.
- 35 Stergiou GS, Karpettas N, Protogerou A, Nasothimiou EG, Kyriakidis M. Diagnostic accuracy of a home blood pressure monitor to detect atrial fibrillation. *J Hum Hypertens*. 2009 Oct;23(10):654-8.
- 36 Kearley K, Selwood M, Van den Bruel A, Thompson M, Mant D, Hobbs FR, Fitzmaurice D, Heneghan C. Triage tests for identifying atrial fibrillation in primary care: a diagnostic accuracy study comparing single-lead ECG and modified BP monitors. *BMJ Open*. 2014 May 2;4(5):e004565.
- 37 Wiesel J, Wiesel DJ, Messineo FC. Home monitoring with a modified automatic sphygmomanometer to detect recurrent atrial fibrillation. *J Stroke Cerebrovasc Dis*. 2007 Jan-Feb;16(1):8-13.
- 38 Wiesel J, Abraham S, Messineo FC. Screening for asymptomatic atrial fibrillation while monitoring the blood pressure at home: trial of regular versus irregular pulse for prevention of stroke (TRIPPS 2.0). *Am J Cardiol*. 2013 Jun 1;111(11):1598-601.
- 39 Ermini, G., et al., Switching from traditional to automatic sphygmomanometer increase opportunistic detection of atrial fibrillation in hypertensive patients. *BJMP*, 2013. 6: p. a6161.
- 40 EREWASH Clinical Commissioning Group, Derbyshire. AF Detection Programme: Quality, Innovation, Prevention & Productivity (QIPP) in Action. 2013. <http://www.atrialfibrillation.org.uk/files/file/Clinicians%20Area/130408-sh-2-Erewash%20case%20study.pdf>
- 41 Kollias A, Xilomenos A, Protogerou A, Dimakakos E, Stergiou GS. Automated determination of the ankle-brachial index using an oscillometric blood pressure monitor: validation vs. Doppler measurement and cardiovascular risk factor profile. *Hypertens Res*. 2011 Jul;34(7):825-30.
- 42 Perk J, De Backer G, Gohlke H, Graham I, Reiner Z, Verschuren M, Albus C, Benlian P, Boysen G, Cifkova R, Deaton C, Ebrahim S, Fisher M, Germano G, Hobbs R, Hoes A, Karadeniz S, Mezzani A, Prescott E, Ryden L, Scherer M, Syv anne M, Scholte op Reimer WJ, Vrints C, Wood D, Zamorano JL, Zannad F; European Association for Cardiovascular Prevention & Rehabilitation (EACPR); ESC Committee for Practice Guidelines (CPG). European Guidelines on cardiovascular disease prevention in clinical practice (version 2012). The Fifth Joint Task Force of the European Society of Cardiology and Other Societies on Cardiovascular Disease Prevention in Clinical Practice (constituted by representatives of nine societies and by invited experts). *Eur Heart J*. 2012 Jul;33(13):1635-701.
- 43 Nasothimiou E, Kollias A, Dafni A, Tzamouranis D, Roussias GS. Performance of 24-hour ambulatory blood pressure monitoring in atrial fibrillation *J Hypertens* 2013;31(e-Supplement A):e495.
- 44 Dafni M, Nasothimiou N, Triantafyllidou E, Kalogeropoulos P, Stergiou GS. Automated atrial fibrillation diagnosis during oscillometric 24-hour ambulatory blood pressure monitoring: a pilot controlled study. *J Hypertens* 2013;31(e-Supplement A):e492.
- 45 Kollias A, Dafni M, Kalogeropoulos P, Nasothimiou E, Destounis A, Tzamouranis D, Roussias L, Stergiou GS. 24-hour ambulatory blood pressure variability in patients with atrial fibrillation. *J Hypertens* 2013;31(e-Supplement A):e505.
- 46 Stergiou GS, Nasothimiou EG, Destounis A, Poulidakis E, Evagelou I, Tzamouranis D. Assessment of the diurnal blood pressure profile and detection of non-dippers based on home or ambulatory monitoring. *Am J Hypertens*. 2012 Sep;25(9):974-8.
- 47 Campbell NR, Berbari AE, Cloutier L, Gelfer M, Kenerson JG, Khalsa TK, Lackland DT, Lemogoum D, Mangat BK, Mohan S, Myers MG, Niebylski ML, O'Brien E, Stergiou GS, Veiga EV, Zhang XH. Policy Statement of the World Hypertension League on Noninvasive Blood Pressure Measurement Devices and Blood Pressure Measurement in the Clinical or Community Setting. *J Clin Hypertens (Greenwich)*. 2014 May 8.
- 48 Wiesel J, Arbesfeld B, Schechter D. Comparison of the microlife blood pressure monitor with the omron blood pressure monitor for detecting atrial fibrillation. *Am J Cardiol* 2014;114:1046-8.
- 49 Gandolfo C., Balestrino M., Bruno C., Finocchi C., Reale N., Validation of a simple method for atrial fibrillation screening in patients with stroke *Neuro Sci*. 2015, April 30.
- 50 Chan PH, Wong CK, Pun L, Wong YF, Wong MM, Chu DW, Siu CW. Head-to-Head Comparison of the AliveCor Heart Monitor and Microlife WatchBP Office AFIB for Atrial Fibrillation Screening in a Primary Care Setting. *Circulation*. 2017 Jan 3;135(1):110-112. doi: 10.1161/CIRCULATIONAHA.116.024439. PubMed PMID: 28028066.
- 51 Chan PH, Wong CK, Pun L, Wong YF, Wong MM, Chu DW, Siu CW. Diagnostic performance of an automatic blood pressure measurement device, Microlife WatchBP Home A, for atrial fibrillation screening in a real-world primary care setting. *BMJ Open*. 2017 Jun 15;7(6):e013685. doi: 10.1136/bmjopen-2016-013685. PubMed PMID: 28619766; PubMed Central PMCID: PMC5577883.

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di Planet Pharma S.p.A. con sede in Cinisello Balsamo (MI), Via Aquileia n. 37, - C.F. 04491260156